

## Attachable Trolley 2

The space-saving transport solution

The handling of equipment has to be a well-practised matter of routine during everyday procedures at hospitals and during transports. Everything necessary also has to be within handy reach. Therefore is now the Attachable Trolley 2, which attaches to patient beds in a flexible and space-saving fashion.

- A well thought out system for transport
- Adds only 20 centimetres [less than 8 inches] in height to the equipment
- Quick and simple to attach
- Simple changing of swivelling oxygen cylinders
- Very stable base
- Tall universal rail for lots of accessories
- Modular system: for custom made solutions



## Much equipment - hands free

When every second counts, it is easy to attach the unit to the patient. Connect the trolley easy and tightly on the bed bar or on the head of the bed. This way the bed can be put in a shock position during transports, without problems .



## A spatial miracle with safety

Since the Attachable Trolley 2 only lengthens patient bed by approx. 20cm [less than eight inches] (even when equipped with an Oxylog 2000 or 3000),

there is no longer anything in the way of even transports in elevators. Its high degree of stability against tipping over provides additional safety during transports.



## Flexible and modular

You decide what the Attachable Trolley takes along with it. Space has been created in a space-saving fashion for all necessary equipment. A transport respirator, such as the Oxylog 3000, a transport monitoring system, syringe pumps along with infusion bags and an aspirator system for patient care can all be easily accommodated. There is an easy changing swivel-out oxygen cylinder. It is also possible to take along gas cylinders of various sizes (e.g. one ten-litre or one five-litre or two three-litre cylinders). For this helpful attendant, it is no problem to add more extensions, for example, a tray, baskets or a Dräger six-place multiple-socket outlet. Its long standard rail provides for a high degree of flexibility when selecting accessories.

## Order List

Designation	Order ID
<b>Attachable Trolley 2</b> Comes with:	MP00011
• Trolley, including belt coupling adjustable in height	
• Standard rail adjustable in height for accommodating transport respirator equipment and other rail-mounted accessories	
• A second standard rail adjustable in height for accommodating a bronchial aspirator system	
• Plastic cylinder jacket	
• Cylinder rack for 10-11 litre compressed-air cylinder	
• Infusion bottle rack for four infusion bottles	
<b>Accessories available:</b>	
• Second cylinder rack for 10-11 litre compressed-air cylinder	MP00043
• Adaptor for two 2-litre or 3-litre cylinders (110 mm diameter)	MP00045
• Rail for aspirator system (short rail)	MP00046
• Rail for emergency respiration equipment (long bar)	MP00047
• Infusion bottle bar, for four units	MP00049
• Infusion bottle rack, for four units	MP00050
<b>Additional accessories:</b> Optionally available on request	

Europe, Middle East, Africa, Latin America, Asia, Pacific:

**Dräger Medical AG & Co. KGaA**  
Moislinger Allee 53-55  
23542 Lübeck  
GERMANY  
Tel: +49-1805-3 72 34 37  
+49-451-882-808  
Fax: +49-451-882-37 79  
E-mail:  
Business.Support@draeger.com

USA:

**Draeger Medical, Inc.**  
3135 Quarry Road  
Telford, PA 18969  
USA  
Tel: +1-215-721-5400  
Toll-free: +1-800-437-2437  
Fax: +1-215-723-5935  
E-mail:  
www.info@draegermed.com

Canada:

**Draeger Medical Canada Inc.**  
120 East Beaver Creek Road  
Suite 104  
Richmond Hill Ontario L4B 4V1  
CANADA  
Tel: +1-905-763-3702  
Toll-free: +1-866-343-2273  
Fax: +1-905-763-1890  
E-mail:  
Canada.Support@draeger.com

www.draeger-medical.com

Manufacturer: Dräger Medical b.v., 5683 Best, The Netherlands

The quality management system at Dräger Medical b.v. is certified according to ISO 13485, ISO 9001 and Annex II of Directive 93/42/EEC (Medical devices).